

ID	Height	Age	Gender	Test Date / Time
John Doe Jr.	139.4cm	10	Male	03.31.2025 16:40

Body Composition Analysis

Total amount of water in my body	Total Body Water	(L)	18.9 (18.0 ~ 22.0)
What I need to build muscles	Protein	(kg)	5.0 (4.9 ~ 5.9)
What I need for strong bones	Minerals	(kg)	1.80 (1.66 ~ 2.04)
Where my excess energy is stored	Body Fat Mass	(kg)	9.3 (3.8 ~ 7.7)
Sum of the above	Weight	(kg)	35.0 (27.3 ~ 36.9)


Muscle-Fat Analysis

	Under	Normal	Over
Weight (kg)	55 70 85 100 115 130 145 160 175 190 205 %	35.0	
SMM (kg)	70 80 90 100 110 120 130 140 150 160 170 %	13.0	
Body Fat mass (kg)	40 60 80 100 160 220 280 340 400 460 520 %	9.3	

Obesity Analysis

	Under	Normal	Over
BMI (kg/m ²)	7.9 10.9 13.9 16.4 18.6 20.2 22.2 24.2 26.2 28.2 30.2	18.0	
PBF (%)	0.0 5.0 10.0 15.0 20.0 25.0 35.0 40.0 45.0 50.0	26.5	

Growth Graph



Body Composition History

Height (cm)	134.5	135.2	136.4	137.2	137.9	138.5	139.0	139.4
Weight (kg)	30.8	31.3	32.0	32.8	33.5	34.0	34.4	35.0
SMM (kg)	12.5	12.7	12.8	13.0	13.1	13.1	13.2	13.0
PBF (%)	20.4	20.7	21.6	22.3	23.1	24.3	25.1	26.5
☒ Recent ☐ Total	07.15.23 14:22	11.19.23 09:30	01.29.24 15:18	03.15.24 11:00	06.21.24 15:00	09.19.24 14:52	12.20.24 15:12	03.31.25 16:40

Growth Score

85 / 100 Points

* If tall and within great body comparison standards, the growth score may surpass 100 points.

Nutrition Evaluation

Protein	☒ Normal	☐ Deficient	
Minerals	☒ Normal	☐ Deficient	
Body Fat	☐ Normal	☐ Deficient	☒ Excessive

Obesity Evaluation

BMI	☒ Normal	☐ Under	☐ Slightly Over	☐ Over
PBF	☐ Normal	☐ Slightly Over	☒ Over	

Body Balance Evaluation

Upper	☒ Balanced	☐ Slightly Unbalanced	☐ Extremely Unbalanced
Lower	☒ Balanced	☐ Slightly Unbalanced	☐ Extremely Unbalanced
Upper-Lower	☒ Balanced	☐ Slightly Unbalanced	☐ Extremely Unbalanced

Segmental Lean Analysis

Right Arm	0.95 kg
Left Arm	0.94 kg
Trunk	10.8 kg
Right Leg	3.41 kg
Left Leg	3.37 kg

Research Parameters

Basal Metabolic Rate	925 kcal (948 ~ 1077)
Child Obesity Degree	109 % (90 ~ 110)

Whole Body Phase Angle

$\phi(^{\circ})$ 50 kHz | 4.3 $^{\circ}$

Segmental Phase Angle

	RA	LA	TR	RL	LL
$\phi(^{\circ})$ 5 kHz	1.4	1.4	3.0	1.9	1.8
50 kHz	3.6	3.3	6.8	5.0	4.8
250 kHz	3.7	3.6	9.4	5.0	4.9

Impedance

